

About this resource

With collaborators & support from multiple institutions, including:



Open Critical Care



USAID
FROM THE AMERICAN PEOPLE



Center for Health Equity
in Surgery & Anesthesia
chsa.ucsf.edu

Last updated April 19, 2021

Disclaimer

This document is intended to be educational in nature and is not a substitute for clinical decision making based on the medical condition presented.

This content is a collaborative effort by representatives of multiple institutions, and this document and the information herein are intended and designed for educational purposes only. You should not rely on this information to replace professional medical advice, diagnosis, and/or treatment, nor should this information ever be used as a substitute for either manufacturers' instructions and guidance or clinical decision-making based on the medical condition presented. It is the responsibility of the user to ensure that all information contained herein is current and accurate, and the creators and hosts of this content make no claims or warranties as to the currency, accuracy, or suitability of this information for any purpose. Any reference to specific equipment, pharmaceuticals, or other medical devices in this document is not meant as an endorsement of such items, and you should consult manufacturers' documentation prior using any such items that may be referenced here. The use of any information in this document is undertaken solely at your own risk, and the creators and hosts of this content shall not be liable for any damages, losses, or other injury caused by the use of any information in this document, nor for any reliance on the accuracy or reliability of such information.

How to Use This Document

This is a living document, created by created by nurses, physicians, respiratory therapists and other healthcare providers from multiple institutions and multiple countries via the OpenCriticalCare.org project.

The goal of this document is to provide tools that can be locally modified to help healthcare providers learning to provide respiratory care for hospitalized patients.

Please **copy or download this file and then modify the document** to fit your needs, the most current data and local resources. For Google documents, this can be done by going to the "File Menu" at the top left of this window and selecting "File" → "Download" then select the format you prefer. Alternatively you can select "File" → "Make a Copy" then use a copy in Google Slides to make your own edits.

To Print: select "Print" from the file menu at the top left. You may need to scale the size to fit your preferred paper size.

Please check back regularly for updates and send any questions or comments to us [here](#).

This work is licensed under a Creative Commons Attribution-NonCommercial-Sharealike 4.0 International License.




Ventilator Circuit Setup

Filter Placement & Humidification Types

“SINGLE LIMB” CIRCUITS*


Heated Wire

(Active heat & humidification)



Heat Moisture Exchanger Filter (HMEF)


HMEFs must always be placed between the circuit wye and the endotracheal tube. They should not be used at the ventilator or exhalation valve.



Heat Moisture Exchanger (HME)

(BV filter by Exhalation Valve)


When using an HME, the expiratory BV filter should ideally be placed either at the exhalation valve (as shown) or at the takeoff of the expiratory limb from the Wye (not shown). This illustration shows a slightly different type of expiratory valve than the others on this page. This has no implication for the location of HMEs or filters.



Heat Moisture Exchanger

(BVF by HME)

When unable to place the expiratory BV filter in the expiratory limb due to hoses glued to the Wye or expiratory valve, it may be placed at the wye in series with the HME. **CAUTION:** this increases dead space & resistance. Review manufacturer's reported dead-space & resistance as tidal volume or pressure support level may require adjustment to ensure adequate ventilation.




* The term 'single limb' is a commonly used misnomer to describe the circuits shown here. These circuits do have a second (expiratory) limb, and thus may be considered functionally similar to 'dual limb circuits'


DUAL LIMB CIRCUITS

Heated Wire


(Active heat & humidification)



Heat Moisture Exchanger (HME)




Heat Moisture Exchange Filter (HMEF)



END TIDAL CO₂ PLACEMENT


Sidestream CO₂

Sampling HMEF with integrated sampling port




Sidestream CO₂

Sampling HME and Bacterial Viral Filter



Mainstream CO₂

Monitoring HME



Deadspace and resistance are increased by the addition of any filters or other adapters, including end tidal CO₂. Read manufacturers' specifications to quantify potential impact on ventilator strategy.


CPAP/NIPPV Circuit Setup

Mask types & filter placement

Always check manufacturer's recommendations to determine if your device, circuit, and mask combination requires an additional anti-asphyxia valve (not shown)

Dual Limb Circuit + Non-Vented Mask


(Exhalation valve in the circuit)



- Because this setup has an **expiratory valve** at or in the ventilator (not shown), a **non-vented mask** must be used for CPAP/NIPPV
- See Figure "Ventilator Circuit Setup" for alternative humidification & filter setups. An HMEF is shown here but may not be the optimal setup depending on local resources.
- Of note, the addition of filters between the circuit Wye and the patient introduce dead space.

"Single Limb" Circuit + Non-Vented Mask


(Exhalation valve in the circuit)



- The term 'single limb' circuit is a commonly used misnomer when referring to the type of circuit shown here. The setup shown has a short expiratory limb and expiratory valve and functions similarly to a dual limb circuit.
- Because this setup has an **expiratory valve**, a **non-vented mask** must be used for CPAP/NIPPV
- See Figure "Ventilator Circuit Setup" for alternative humidification & filter setups.

True Single Limb Circuit + Vented Mask


(No exhalation valve)



- The 'true' single limb circuit is commonly found in home CPAP/NIPPV machines and has **no expiratory valve**.
- Exhaled breath must leave through intentional leak ports - in this setup a **vented mask** is used
- The BV filter in this scenario protects the patient from a potentially contaminated device/circuit but does not protect personnel from potentially infectious aerosolized particles in the exhaled gas.


True Single Limb Circuit + Leak Adapter + Non-Vented Mask

(No exhalation valve)




- This 'true' single limb circuit is commonly found in home CPAP/NIPPV machines and has **no expiratory valve**.
- Exhaled breath must leave through intentional leak ports - in this setup, a non-vented mask is used and the leak occurs via a specialized **fixed leak adapter** in the circuit.
- Either a bacterial viral filter (if active heated humidification is being used) or HMEF (if no active heated humidification) should be placed between the leak adapter & the patient.

Bag valve resuscitator setup



HMEF may be considered over a bacterial viral filter if duration of use on the bag valve resuscitator is anticipated to be prolonged.

In-line suction catheter setup



This illustration shows the location of the in-line suction catheter in relation to the patient circuit. There are other potential configurations not shown here, but in all, the in-line suction catheter is placed on the endotracheal tube.